



Inbios International, Inc.
562 1st Ave. South
Suite 600
Seattle, WA 98104
206-344-5821 (phone)
206-344-5823 (fax)
www.inbios.com

Product Specification Sheet

Product	Monoclonal antibody specific to <i>Campylobacter jejuni</i>
Catalog Number	INB9944
Description	Antibody specific to <i>Campylobacter jejuni</i>
Specificity	
Clone	CCJ1
Host Animal	Mouse
Format	Purified, liquid
Isotype	IgG
Buffer	PBS pH 7.4
Application	Suitable for use in ELISA. Each laboratory should determine an optimum working titer for their own application. Other formats have not been tested.
Storage	Store at -20°C. Avoid multiple freeze/thaw cycles. Prepare working dilution only prior to immediate use.

For Research Use Only. Not for use in Diagnostic Procedures.

Important Note: Centrifuge before opening to ensure complete recovery of vial contents.